## MODUL PRAKTIKUM KIMIA FISIKA FARMASI

Pas foto terkini

**Pemilik buku**

Nama : \_\_\_\_\_  
(huruf cetak)

NRP : \_\_\_\_\_

Kelas : \_\_\_\_\_

Kelompok : \_\_\_\_\_ No. : \_\_\_\_\_

**LABORATORIUM KIMIA FISIKA FARMASI**  
**FAKULTAS FARMASI**  
**UNIVERSITAS SURABAYA**

## KATA PENGANTAR

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4. PUSTAKA

Sinko PJ, 2012, Physical Pharmacy and Pharmaceutical Sciences 6<sup>th</sup> ed., Lippincott Williams and Wilkins.

Remington The Science and Practice of Pharmacy 2011.

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
### Ostwald viscometers

#### Ostwald viscometers

- filling quantity: 3 ml
- overall length: approx. 220 mm

with ring marks,  
without constant,  
for manual measurements

Type No.	Order No.	Capillary Ø i [mm]	Transit time for water approx. [s]	Constant K (approx.)	for use from [mm <sup>2</sup> /s] (approx.)
509 03	285404006	0.3	250	0.004	0.3
509 04	285404014	0.4	75	0.01	1
509 05	285404022	0.5	30	0.03	2.5
509 06	285404039	0.6	15	0.07	5.5
509 07	285404047	0.7	10	0.1	10

## PHARMACEUTICAL SOLVENTS

The remarkable growth of the solvent industry is attested by the more than 300 solvents now being produced on an industrial scale. Chemically, these include a great variety of organic compounds, ranging from hydrocarbons through alcohols, esters, ethers, and acids to nitroparaffins. Their main applications are in industry and the synthesis of organic chemicals. Comparatively few, however, are used as solvents in pharmacy, because of their toxicity, volatility, instability, and/or flammability. Those commonly used as pharmaceutical solvents are described in this section.

### ACETONE

#### 2-Propanone; Dimethyl Ketone

$\text{CH}_3\text{COCH}_3$

Acetone [67-64-1]  $\text{C}_3\text{H}_6\text{O}$  (58.08).

**Caution**—It is very flammable. Do not use where it may be ignited.

**Preparation**—Formerly obtained exclusively from the destructive distillation of wood. The distillate, consisting principally of methanol, acetic acid, and acetone was neutralized with lime, and the acetone was separated from the methyl alcohol by fractional distillation. Additional quantities were obtained by pyrolysis of the calcium acetate formed in the neutralization of the distillate.

It now is obtained largely as a by-product of the butyl alcohol industry. This alcohol is formed in the fermentation of carbohydrates such as corn starch, molasses, etc, by the action of the bacterium *Clostridium acetobutylicum* (Weizmann fermentation), and it is always one of the products formed in the process. It also is obtained by the catalytic oxidation of isopropyl alcohol, which is prepared from propylene resulting from the cracking of crude petroleum.

**Description**—Transparent, colorless, mobile, volatile, flammable liquid with a characteristic odor; specific gravity not more than 0.789;

distills between 55.5 and 57°C; congeals about -95°C; aqueous solution neutral to litmus.

**Solubility**—Miscible with water, alcohol, ether, chloroform, or most volatile oils.

**Uses**—An antiseptic in concentrations above 80%. In combination with alcohol it is used as an antiseptic cleansing solution. It is employed as a menstruum in the preparation of oleoresins in place of ether. It is used as a solvent for dissolving fatty bodies, resins, pyroxylin, mercurials, etc, and also in the manufacture of many organic compounds such as chloroform, chlorobutanol, and ascorbic acid.


### ALCOHOL

#### Ethanol; Spiritus Vini Rectificatus; S. V. R.; Spirit of Wine; Methylcarbinol

Ethyl alcohol [64-17-5]; contains 92.3 to 93.8%, by weight (94.9 to 96.0%, by volume), at 15.56°C (60°F) of  $\text{C}_2\text{H}_5\text{OH}$  (46.07).

**Preparation**—Has been made for centuries by fermentation of certain carbohydrates in the presence of zymase, an enzyme present in yeast cells. Usable carbohydrate-containing materials include molasses, sugar cane, fruit juices, corn, barley, wheat, potato, wood, and waste sulfite liquors. As yeast is capable of fermenting only D-glucose, D-fructose, D-mannose, and D-galactose, it is essential that more complex carbohydrates, such as starch, be converted to one or more of these simple sugars before they can be fermented. This is accomplished variously, commonly by enzyme- or acid-catalyzed hydrolysis.

The net reaction that occurs when a hexose, glucose for example, is fermented to alcohol may be represented as

but the mechanism of the process is very complex. The fermented liquid, containing about 15% alcohol, is distilled to obtain a distillate contain-

ing 94.9% C<sub>2</sub>H<sub>5</sub>OH, by volume. To produce absolute alcohol, various processes dehydrate the 96% product.

Hydration of ethylene, abundant supplies of which are available from natural and coke oven gases, from waste gases of the petroleum industry, and other sources may produce it also. In another synthesis acetylene is hydrated catalytically to acetaldehyde, which then is hydrogenated catalytically to ethyl alcohol.

**Description**—Transparent, colorless, mobile, volatile liquid; slight but characteristic odor, burning taste; boils at 78°C but volatilizes even at a low temperature, and is flammable; when pure, it is neutral toward all indicators; specific gravity at 15.56 (the US Government standard temperature for Alcohol) not above 0.816, indicating not less than 92.9% of C<sub>2</sub>H<sub>5</sub>OH by weight, or 94.9% by volume.

**Solubility**—Miscible with water, acetone, chloroform, ether, or many other organic solvents.

**Incompatibilities**—This and preparations containing a high percentage of alcohol will precipitate many inorganic salts from an aqueous solution. Acacia generally is precipitated from a hydroalcoholic medium when the alcohol content is greater than about 35%.

Strong oxidizing agents such as chlorine, nitric acid, permanganate, or chromate in acid solution react, in some cases violently, with it to produce oxidation products.

Alkali cause a darkening in color because of the small amount of aldehyde usually present in it.

**Uses**—In pharmacy principally for its solvent powers. It also is used as the starting point in the manufacture of many important compounds, like ether, chloroform, etc. It also is used as a fuel, chiefly in the denatured form.

It is a CNS depressant. Consequently, it occasionally has been administered intravenously for preoperative and postoperative sedation for patients in whom other measures are ineffective or contraindicated. The dose employed is 1 to 1.5 mL/kg. Its intravenous use is a specialized procedure and should be employed only by one experienced in the technique of such use.

It is used widely and abused by lay persons as a sedative. It has, however, no medically approved use for this purpose. Moreover, alcohol potentiates the CNS effects of numerous sedative and depressant drugs. Hence, patients taking certain prescription drugs or OTC medications should not use it.

Externally, it has a number of medical uses. It is a solvent for the toxicodendrol causing by poisoning and should be used to wash the skin thoroughly soon after contact. In a concentration of 25% it is employed for bathing the skin for the purpose of cooling and reducing fevers. In high concentrations it is a rubefacient and an ingredient of many liniments. In a concentration of 50% it is used to prevent sweating in astringent and anhidrotic lotions. It also is employed to cleanse and harden the skin and is helpful in preventing bedsores in bedridden patients. In a concentration of 60 to 90% it is germicidal. At optimum concentration (70% by weight) it is a good antiseptic for the skin (local antiseptic) and also for instruments. It also is used as a solvent to cleanse the skin splashed with phenol. High concentrations of it often are injected into nerves and ganglia for the relief of pain, accomplishing this by causing nerve degeneration.

#### DENATURED ALCOHOL

An act of Congress, June 7, 1906, authorizes the withdrawal of alcohol from bond without the payment of internal revenue tax, for the purpose of denaturation and use in the arts and industries. This is ethyl alcohol to which has been added such denaturing materials as to render the alcohol unfit for use as an intoxicating beverage. It is divided into two classes, namely, completely denatured alcohol and specially denatured alcohol, prepared in accordance with approved formulas prescribed in Federal Industrial Alcohol Regulations 3.

**Completely Denatured Alcohol**—This term applies to ethyl alcohol to which has been added materials (methyl isobutyl ketone, pyrene, gasoline, acetone, kerosene, etc.) of such nature that the products may be sold and used within certain limitations without permit and bond.

**Specially Denatured Alcohol**—This alcohol is intended for use in a greater number of specified arts and industries than completely denatured alcohol, and the character of the denaturant or denaturants used is such that specially denatured alcohol may be sold, possessed, and used only by those persons or firms that hold basic permits and are covered by bond.

**Uses**—Approximately 50 specially denatured alcohol formulas containing combinations of more than 90 different denaturants are available to fill the needs of qualified users. Large amounts of specially denatured alcohols are used as raw materials in the production of acetaldehyde, synthetic rubber, vinegar, and ethyl chloride as well as in the manufacture of proprietary solvents and cleaning solutions. Ether and chloroform can be made from suitably denatured alcohols, and fur-

malin for the manufacture of Iodine Tincture, Green Soap Tincture, and Rubbing Alcohol are set forth in the regulations.

Specially denatured alcohols also are used as solvents for surface coatings, plastics, inks, toilet preparations, and external pharmaceuticals. Large quantities are used in the processing of such food and drug products as pectin, vitamins, hormones, antibiotics, alkaloids, and blood products. Other uses include supplemental motor fuel, rocket and jet fuel, antifreeze solutions, refrigerants, and cutting oils. Few products are manufactured today that do not require the use of alcohol at some stage of production. Specially denatured alcohol may not be used in the manufacture of foods or internal medicines when any of the alcohol remains in the finished product.

#### DILUTED ALCOHOL

##### Diluted Ethanol

A mixture of alcohol and water containing 41.0 to 42.0%, by weight (48.4 to 49.5%, by volume), at 15.56°C, of C<sub>2</sub>H<sub>5</sub>OH (46.07).

##### Preparation

Alcohol 500 mL  
Purified Water 500 mL

Measure the alcohol and the purified water separately at the same temperature, and mix. If the water and the alcohol and the resulting mixture are measured at 25°C, the volume of the mixture will be about 970 mL.

When equal volumes of alcohol and water are mixed together, a rise in temperature and a contraction of about 3% in volume take place. In small operations the contraction generally is disregarded; in larger operations it is very important. If 50 gal of official alcohol are mixed with 50 gal of water, the product will not be 100 gal of diluted alcohol, but only 96 1/4 gal, a contraction of 3 3/4 gal. US Proof Spirit differs from this and is stronger; it contains 50%, by volume, of absolute alcohol at 15.56°C (60°F). This corresponds to 42.5% by weight and has a specific gravity of 0.9341 at the same temperature. If spirits have a specific gravity lower than that of proof spirit (0.9341), they are said to be above proof; if greater, below proof.

It also may be prepared from the following:

Alcohol 408 g  
Purified Water 500 g

**Rules for Dilution**—The following rules are applied when making an alcohol of any required lower percentage from an alcohol of any given higher percentage:

**I. By Volume**—Designate the volume percentage of the stronger alcohol by *V* and that of the weaker alcohol by *v*.

**Rule**—Mix *v* volumes of the stronger alcohol with purified water to make *V* volumes of product. Allow the mixture to stand until full contraction has taken place and until it has cooled, then make up the deficiency in the *V* volumes by adding more purified water.

**Example**—An alcohol of 30% by volume is to be made from an alcohol of 94.9% by volume.—Take 30 volumes of the 94.9% alcohol, and add enough purified water to produce 94.9 volumes at room temperature.

**II. By Weight**—Designate the weight-percentage of the stronger alcohol by *W* and that of the weaker alcohol by *w*.

**Rule**—Mix *w* parts by weight of the stronger alcohol with purified water to make *W* parts by weight of product.

**Example**—An alcohol of 50% by weight is to be made from an alcohol of 92.3% by weight.—Take 50 parts by weight of the 92.3% alcohol, and add enough purified water to produce 92.3 parts by weight.

**Description**—As for Alcohol except its specific gravity is 0.935 to 0.937 at 15.56°C, indicating that the strength of C<sub>2</sub>H<sub>5</sub>OH corresponds to that given in the official definition.

**Uses**—A solvent in making tinctures, fluid-extracts, extracts, etc. Its properties already have been described fully in connection with the various preparations. Its value consists not only in its antiseptic properties, but also in its possessing the solvent powers of both water and alcohol. See Alcohol.

**CHLOROFORM**—page 1085.

#### GLYCERIN

##### 1,2,3-Propanetriol Glycerol

Glycerol [56-81-5] C<sub>3</sub>H<sub>8</sub>O<sub>3</sub> (92.09).

Chemically, it is the simplest trihydric alcohol. It is worthy of special note because the two terminal alcohol groups are primary, whereas the middle one is secondary. Thus this becomes the first polyhydric alcohol that can yield both an aldose (glyceraldehyde) and a ketose (dihydroxyacetone).

# Viskometer Bola Jatuh (*Falling Sphere Viscometer*)

## Unit Description

---

### 4. Unit Description

#### 4.1 Principle of the measurement

The HAAKE Falling Ball Viscometer measures the viscosity of transparent newtonian liquids. This viscosity is correlated to the time a ball requires to fall a defined distance. The rolling and sliding movement of the ball through the sample filled into a slightly inclined cylindrical measuring tube is described by means of the fall time. The test results are given as the dynamic viscosity using the internationally standardized absolute unit of "milli Pascal seconds" (mPa·s).


Note: 1 mPa·s = 1 cP (centi Poise)

The Falling Ball Viscometer corresponds to the requirements of many international standards, i.e. ISO 12058 and the german standard DIN 53 015.

#### 4.2 Description of the instrument

The heart of the instrument is the measuring tube made of glass **1** and a ball **2**. This tube carries two ring marks A and B, which are spaced 100 mm apart and which limit the measuring distance (ring mark C is equidistant between A and B). The measuring tube is jacketed by means of an outer glass tube, which encloses a room **3** to be filled with a temperature controlled liquid. The measuring tube is fastened to the stand in such a way that its axis is inclined with respect to the vertical by 10° during the measurement.


The measuring tube together with the jacket may be pivoted in order to turn the tube upside down again to let the ball return to the initial position before a measurement. The measuring tube is closed on both sides by two stoppers, one of which **13** contains a capillary and a small reservoir. This stopper prevents undesirable changes of pressure in the liquid sample and has a passage for air bubbles when the temperature is being changed. The viscometer incloses all samples completely to prevent volatization and film forming. The stand may be levelled by means of its water level **16** and the levelling screws **17**. The easily interchangeable thermometer allows a precise temperature control.

## Functional Elements

---

### 5. Functional Elements

- |                  |                            |
|------------------|----------------------------|
| 1 Falling tube   | 12 Jacket tube             |
| 2 Ball           | 13 Hollow stopper          |
| 3 Tempering room | 14 Stopper                 |
| 4 Screw          | 15 Capillary               |
| 7 Set screw      | 16 Closing plate           |
| 8 Gasket         | 17 Gasket for falling tube |
| 9 Cover          | 18 Threaded bush           |
| 11 Brace         | 20 Connecting rod          |

# Measuring

---

## 6. Measuring

### 6.1 Preparation for a test

#### 6.1.1 Temperature control

The Falling Ball Viscometer may be temperature controlled in a temperature range from  $-20$  up to  $+120^{\circ}\text{C}$  using liquid circulators, i.e. one of the HAAKE "DC" or "F" series.

The sample should rest at least some 15 minutes in the measuring tube at the test temperature before the measurement is started.

The temperature in the jacket around the measuring tube must be maintained within a temperature tolerance of  $\pm 0.03^{\circ}\text{C}$  for test temperatures between  $10$  up to  $80^{\circ}\text{C}$ . For test temperatures beyond these limits the tolerances may be increased to  $+0.05^{\circ}\text{C}$ .


The tempering room 3 must be free of air bubbles.

#### 6.1.2 Loading the sample

All parts of the viscometer being in direct contact with the sample must be kept clean and dry.

A sample volume of approximately  $45\text{ cm}^3$  is poured into the measuring tube 1 up to 20 mm below the rim of the tube. Then the ball 2 is placed into the tube and the hollow stopper 13 is introduced. The liquid should reach a level just beyond the capillary 15. The sample in the tube must be free of air bubbles.

Before the final test data are taken the ball should run through the tube up and down at least once to improve the homogeneity of the samples and its temperature uniformity.

### 6.1.3 Selection of the balls

The standard ball set contains 6 balls, which pass through the measuring tube of an inner diameter of approximately  $15.94 \pm 0.01$  mm.

Order-No.	Ball No.	Made of	Density $\rho$ g/cm <sup>3</sup>	Diameter of the ball mm	Constant K (approx.) mPa·s·cm <sup>3</sup> /g·s	Recomm. measuring range mPa·s
800-0002	1	boron silica glass	2.2	$15.81 \pm 0.01$	0.007	0.6 – 10
800-0003	2	boron silica glass	2.2	$15.6 \pm 0.05$	0.09	7 – 130
800-0004	3	nickel iron alloy	8.1	$15.6 \pm 0.05$	0.09	30 – 700
800-0005	4	nickel iron alloy	8.1	$15.2 \pm 0.1$	0.7	200 – 4800
800-0006	5	W.-No. 4034	7.7–8.1	$14.0 \pm 0.5$	4.5	800 – 10000
800-0007	6	W.-No. 4034	7.7–8.1	$11.0 \pm 1$	33	6000 – 75000

Additionally the following ball is deliverable:

Order-No.	Ball	Made of	Density $\rho$ g/cm <sup>3</sup>	Diameter of the ball mm	Constant K (approx.) mPa·s·cm <sup>3</sup> /g·s	Recomm. measuring range mPa·s
800-0009	G	boron silica glass	2.2	$15.91 \pm 0.01$	–	gases

The measuring ranges for viscosity indicated are related to DIN 53015 / ISO 12058.

Sometimes it may be necessary to use two different balls in order to cover a wider measuring range, i.e. when the function of viscosity versus temperature is measured over a wide temperature intervall; in this case, you insert two different balls at the same time into the measuring tube with the smaller ball inserted first. The reduced starting distance will not influence the test result significantly, but increases the uncertainty.

## ***Measuring***

---

### **6.2 Measurement of the falling times**

The jacket tube snaps into a defined 10°-position at the bottom of the instrument.

By turning over the jacket tube, the ball is set to the measuring position.

The falling time of the ball moving from the ring mark A to ring mark B is determined by using a stop watch. The time period starts when the lower periphery of the ball touches the ring mark A, which must appear as a straight line. The falling time ends when the lower periphery of the ball touches the ring mark B, which again must appear as a straight line. If one uses the distance AC or CB to reduce very long falling times for high viscous liquids the double of the measuring time period must be taken into account.

Turning the jacket tube 180° again the ball returns to its start position. It is good practice to take the mean value out of several falling time values (3 to 5).

The falling times for the ball returning may vary from the normal value (up to 1 %). If the returning of the ball should also be used for exact measurements a new constant K must be determined.

**Constant for the returning of the ball:**

$$K_{\text{return}} = \frac{\text{normal falling time} \cdot \text{normal constant } K}{\text{falling time when returning}}$$

When testing dark liquids it is usually very difficult to see the lower part of the ball. In this case we advise to take the ball equator when it passes through the ring marks.

## **Measuring**

---

### **6.3 Evaluation of the test results**

The dynamic viscosity  $\eta$  (in mPa·s) is calculated using the following equation:

$$\eta = K (\rho_1 - \rho_2) \cdot t$$

where:

$K$  = ball constant in mPa·s·cm<sup>3</sup>/g·s (see chapter 11)

$\rho_1$  = density of the ball in g/cm<sup>3</sup> (see chapter 11)

$\rho_2$  = density of the liquid to be measured at the measuring temperature in g/cm<sup>3</sup>

$t$  = falling time of the ball in seconds.

#### **Test results:**

The dynamic viscosity  $\eta$  is given in units of mPa·s (cP) and must be completed by stating the sample temperature.

The dynamic viscosity  $\eta$  may be converted to the kinematic viscosity  $\nu$  by using the following equation:

$$\nu = \frac{\eta}{\rho}$$

$\nu$  = kinematic viscosity [mm<sup>2</sup>/s] [1 mm<sup>2</sup>/s = 1 cSt]

$\eta$  = dynamic viscosity [mPa·s]

$\rho$  = density of the liquid sample [g/cm<sup>3</sup>]

# Viskometer *Cup and Bob (Stormer)*

## Experiment 1

# Viscosity Measurement

### Purpose

The purpose of this experiment is to measure the viscosity of a glycerin-water mixture with a Thomas-Stormer viscometer.

### Apparatus

Figure 1.1 is a schematic of the viscometer. A weight,  $W$ , is used to drive a rotor that is partially submerged in a sample of liquid. The torque exerted by viscous shear on the rotor is balanced by the work input of the falling weight. The experiment involves measurement of the time it takes for a known number of revolutions of the rotor.


Figure 1.1: Thomas-Stormer viscometer.

## Theory

The operation of the viscometer relies on a linear velocity profile in the gap between the rotor and the fixed cylinder. If the velocity profile is linear the viscous shear stress on the surface of the rotor can be written

$$\tau = \mu k \omega \quad (1.1)$$

where  $\tau$  is the viscous shear stress,  $\mu$  is the fluid viscosity,  $k$  is a constant that depends only on the geometry of the viscometer, and  $\omega$  is the angular velocity of the rotor. Given the shear stress from Equation (1.1), the torque exerted by the rotor on the fluid is

$$T_f = (\tau A) r_r \quad (1.2)$$

where  $A$  is the wetted surface area of the rotor, and  $r_r$  is the radius of the rotor. The area,  $A$ , accounts for the inner and outer surfaces of the rotor. Since the fluid is in contact with both surfaces  $r_r$  is an effective radius.

Neglecting any friction in the pulleys and bearings, the power dissipated by viscous stresses in the fluid,  $P_f$ , is equal to the power input of the falling weight,  $P_w$ .

$$P_f = P_w \implies T_f \omega = W V \quad (1.3)$$

where  $W$  is the magnitude of the weight, and  $V$  is the velocity of the falling weight. Combining Equations (1.1) through (1.3) gives

$$\mu k \omega^2 A r_r = W V \quad (1.4)$$

The velocity of the weight falling a distance  $L$  in time  $t$  is

$$V = \frac{L}{t} = \frac{2\pi r_s n_s}{t} \quad (1.5)$$

where  $n_s$  is the number of revolutions of the spindle, and  $r_s$  is the radius of the spindle. The angular velocity of the rotor is

$$\omega = \frac{2\pi n_r}{t} \quad (1.6)$$

where  $n_r$  is the number of revolutions of the rotor in time  $t$ .

Substitution of Equations (1.5) and (1.6) into Equation (1.4) and rearranging yields

$$\mu = \frac{1}{k \omega \frac{2\pi n_r}{t} A r_r} W \frac{2\pi r_s n_s}{t} = \frac{r_s \beta}{k A r_r} \frac{W}{\omega} \quad (1.7)$$

where  $\beta = n_s/n_r$  is the overall gear ratio between the spindle and the rotor. Defining the viscometer constant as

$$C = \frac{r_s \beta}{k A r_r} \quad (1.8)$$

Equation (1.7) can be written

$$\mu = C \frac{W}{\omega} \quad (1.9)$$

Assuming that the model of viscous shear in Equation (1.1) is valid, the constant  $C$  depends only on the geometry of the device.

## Procedure

Use of the Thomas-Stormer viscometer requires determination of  $C$  in Equation (1.9) by calibrating the instrument with a fluid having a known viscosity. With  $C$  known from the calibration step, measurements of  $W$  and  $\omega$  can be used with Equation (1.9) to compute the viscosity of an unknown fluid. The overall procedure may be divided into three phases: (1) setting up the viscometer, (2) adjusting the weight in preparation for the tests, and (3) running the tests.

### Set Up the Viscometer

1. Fill the test cup to the top of the side vanes. Make sure that the depth of the fluid is the same for all tests.
2. Replace the test cup in the viscometer
3. Raise the platform that supports the bath and test cup until the fluid is at least 0.6 cm (0.25 inch) above the top of the rotor. Add fluid if necessary. Make sure that the platform is raised so that it touches the stop.
4. Place the thermometer in the holder. Allow the system to come into thermal equilibrium and record the temperature. Note that most samples have been sitting overnight in the lab. Unless the fluid temperature has been disturbed, e.g., by excessive handling, it should be in reasonable thermal equilibrium with the viscometer and the surroundings.

### Adjust Weight in Preparation for Tests

The viscosity measurement is based on the assumption that the flow on the surface of the rotor is laminar. After placing a new sample of liquid in the test cup, and raising the cup into position, adjust the weight until 100 revolutions on the counter takes at least 20 seconds. Shorter run times will cause turbulent flow, and result in erroneously high viscosity values. The preliminary weight adjustment will also allow you to become familiar with the measurement procedure.

All data runs should be taken over a range of weights no greater than the weight that gives 100 revolutions in no less than 20 seconds. You will also need to allow about 20 revolutions for the rotor to attain steady state velocity. Make sure, therefore, that the weight can fall far enough to cause at least 120 revolutions of the rotor.

### Running Tests

Once the instrument is set up, the calibration runs and the viscosity measurement runs use the following procedure.

1. Turn the brake on and raise the driving weight by turning the handle of the rewinding drum counter-clockwise until the weight nearly touches the pulley.
2. Release the brake one quarter turn and allow the weight to slowly descend until the pointer on the dial is located between 80 and 90. Reset the brake. The time for the weight to fall 100 revolutions will be measured. The starting time is when the revolution counter crosses 0, and the stopping time is when the revolution counter crosses 100. Setting the starting position so that the revolution counter is between 80 and 90 allows several revolutions of the rotor to be completed before the beginning of the measured time interval. This guarantees that the rotor is rotating at steady angular velocity during the interval of the time measurement.
3. Reset the stopwatch.
4. Fully release the brake. Start the stopwatch when the counter passes 0 and stop it when the counter passes 100. Record the time and mass of the weight box.
5. Repeat the preceding steps for several weights (at least 5) by adding or removing shot from the weight box.

# Rotoviskometer (*Brookfield Dial*)

## I. INTRODUCTION

The Brookfield Dial Viscometer measures fluid viscosity at given shear rates. Viscosity is a measure of a fluid's resistance to flow. You will find a detailed description of the mathematics of viscosity in the Brookfield publication *"More Solutions to Sticky Problems"* a copy of which was included with your Dial Viscometer. The Dial Viscometer rotates a sensing element in a fluid and measures the torque necessary to overcome the viscous resistance to the induced movement. This is accomplished by driving the immersed element, which is called a spindle, through a beryllium copper spring. The degree to which the spring is wound, indicated by the red pointer, is proportional to the viscosity of the fluid.

The Viscometer is able to measure over a number of ranges since, for a given spring deflection, the actual viscosity is proportional to the spindle speed and is related to the spindle's size and shape. For a material of given viscosity, the resistance will be greater as the spindle size and/or rotational speed increase. The minimum viscosity range is obtained by using the largest spindle at the highest speed; the maximum range by using the smallest spindle at the slowest speed.

There are four basic spring torque series offered by Brookfield:

<u>Model</u>	<u>Spring Torque (Dyne-cm)</u>
LV	673.7
RV	7,187.0
HA	14,374.0
HB	57,496.0

The higher the torque calibration, the higher the measurement range. The viscosity measurement range for each torque calibration may be found in Appendix B.

All units of measurement are calculated in units of centipoise (cP) by using a look-up table known as "the Brookfield Factor Finder" to convert torque reading. The equivalent units of measurement in the SI system are calculated using the following conversions:

Viscosity:	$\frac{\text{CGS}}{1 \text{ cP}}$	=	$\frac{\text{SI}}{1 \text{ mPa}\cdot\text{s}}$	mPa·s = milli Pascal - second
Torque:	$10^7 \text{ dyne-cm}$	=	1 Newton-m	

### I.1 Components

- 1) Dial Viscometer
- 2) Model A Laboratory Stand
- 3) Spindle Set with Case (4 spindles for LVF, LVT; 7 Spindles for RVF, RVT, HAT, and HBT).  
For Cone/Plate versions: a spindle wrench, one cone spindle and a sample cup (Part No. CP-44Y) replace the spindle set
- 4) Guard Leg (LVF, LVT, RVF, and RVT only)
- 5) Carrying Case
- 6) Operating Instruction Manual
- 7) Factor Finder

Please check to be sure that you have received all components, and that there is no damage. If you are missing any parts, please notify Brookfield Engineering or your local Brookfield agent immediately. Any shipping damage must be reported to the carrier.

## 1.2 Utilities

Input Voltage:	115 VAC or 230 VAC
Input Frequency:	50 or 60 Hz
Power Consumption:	Less than 20 WATTS

## 1.3 Specifications

Speeds:	LVF:	60, 30, 12, 6	
	LVT:	60, 30, 12, 6, 3, 1.5, 0.6, 0.3	
	RVF:	20, 10, 4, 2	
	RVT:	100, 50, 20, 10, 5, 2.5, 1, 0.5	
	HAT:	100, 50, 20, 10, 5, 2.5, 1, 0.5	
	HBT:	100, 50, 20, 10, 5, 2.5, 1, 0.5	
Accuracy:	± 1% of Full Scale Range in use (See Appendix D for details)		
Repeatability:	± 0.2% of Full Scale Range in use		
Weight:	Gross Weight	20 lb	9 kg
	Net Weight	17 lb	7.7 kg
	Carton Volume	1.65 cu ft	0.05 m <sup>3</sup>
Operating Temperature:	-20°C (-4°F) to 65°C (149°F)		

## 1.4 Set-Up

- 1) Assemble the laboratory stand. (Refer to Appendix E.)
- 2) Mount the Viscometer securely on a Brookfield laboratory stand. On some Viscometers, it may be necessary to unscrew the nut located at the point where the power cord enters the viscometer. This permits the metal handle to be inserted into the laboratory stand clamp.

**Note:** The position of the laboratory stand clamp assembly is important. Refer to Appendix E for proper alignment and positioning of the clamp assembly.

Level the viscometer, referring to the bubble level on the instrument. If the viscometer cannot be leveled, recheck the laboratory stand assembly, as shown on Sheet #82-0330 (supplied with Model A laboratory stand).

- 3) Verify that the viscometer's power requirements match your power source before connecting it to power.
- 4) For Cone/Plate models, refer to Appendix A.

## II. GETTING STARTED

### II.1 Operation

The Brookfield Viscometer is powered by a precision synchronous motor. Exact speeds of rotation are assured as the motor will turn erratically and spasmodically if synchronization cannot be maintained.

Speed changes are affected by a gear train having either four or eight speeds. Maximum speed (RPM) will be set at full clockwise rotation and minimum speed at full counter-clockwise rotation. Although not absolutely necessary, it is advisable to change speeds while the motor is running.

Four speed Viscometers have a square speed control knob with one number shown on each of four faces. The instrument's rotational speed is indicated by the number facing up. Eight speed models have a square speed control knob with two numbers on each face; by moving the knob through two complete turns, speeds may be changed in sequence. No trouble will be experienced in differentiating between the two speeds shown on each face since each pair is in the ratio of 20:1. To insure rotation at the indicated speed, it is important that the face of the knob upon which this speed is shown be closely parallel to the viscometer's dial.

LV Viscometers are provided with a set of four spindles and a narrow guard leg; RV Viscometers come with a set of seven spindles and a wider guard leg; HA and HB Viscometers come with a set of seven spindles and no guard leg.

The spindles are attached to the viscometer by screwing them to the lower shaft. **Note that the spindles have a left-hand thread.** The lower shaft should be held in one hand and lifted up. The spindle should be screwed to the left. The face of the spindle nut and the matching surface on the lower shaft should be smooth and clean to prevent eccentric rotation of the spindle. Spindles can be identified by the number on the side of the spindle nut.

All Brookfield Dial Reading Viscometers are provided with a clutch lever located at the back of the instrument. Depressing the lever raises the dial against the pointer and "holds" the instrument's reading. When the clutch is released the dial is lowered and the pointer will be freed.

Any of the controls on the Viscometer - the motor switch, speed change knob, and clutch - may be operated independently of the other.

### II.2 Viscometer Readings and Viscosity Calculations

1. Mount the guard leg, if used, (spindle protector) on the viscometer.
2. Attach the spindle (left-hand thread) to the viscometer lower shaft by lifting the coupling screw slightly. Hold it firmly with one hand while screwing the spindle on with the other (note left-hand thread). Avoid putting side thrust on the shaft.

To avoid trapping air bubbles under the disc-type spindles, do the following before attaching the spindle. Begin by immersing the spindle in a diagonal path, slowly drag the spindle across the

fluid surface, and bring the spindle to an upright position and thread onto screw.

3. Lower and center spindle in the test material (600 ml beaker) until the "meniscus" of the fluid is at the center of the immersion groove on the spindle's shaft.
4. To make a viscosity measurement, turn the motor switch "ON". This energizes the viscometer drive motor. Allow time for the indicated reading to stabilize. The time required for stabilization will depend on the speed at which the viscometer is running and the characteristics of the sample fluid.

When making a viscosity measurement, the reading should be noted and multiplied by the factor appropriate to the viscometer model/spindle/speed combination being used. Refer to the table below or to the FACTOR FINDER for calculating viscosity. For maximum accuracy, readings below 10.0% torque (dial reading) should be avoided.

5. Turn the viscometer motor switch "OFF" when changing a spindle, changing samples, etc. Remove spindle before cleaning.
6. The interpretation of results and the instrument's use with non-Newtonian and thixotropic materials is discussed in the booklet, "More Solutions to Sticky Problems" and in Appendix C.

The following tables apply to Brookfield Viscometer Models LV, RV, HA and HB with standard spindles. They enable the user to convert the percent scale reading into a viscosity value in units of centipoise.

To convert the viscometer dial reading to a viscosity value in units of centipoise, multiply the reading noted on dial viscometer by the appropriate factor in the following tables.

LV Series Viscometer							
1		2		3		4	
.3	200	.3	1K	.3	4K	.3	20K
.6	100	.6	500	.6	2K	.6	10K
1.5	40	1.5	200	1.5	800	1.5	4K
3	20	3	100	3	400	3	2K
6	10	6	50	6	200	6	1K
12	5	12	25	12	100	12	500
30	2	30	10	30	40	30	200
60	1	60	5	60	20	60	100

= Spindle       = Spindle Speed  
 = Factor      K = 1000

RV Series Viscometer													
1		2		3		4		5		6		7	
.5	200	.5	800	.5	2k	.5	4k	.5	8k	.5	20k	.5	80k
1	100	1	400	1	1k	1	2k	1	4k	1	10k	1	40k
2	50	2	200	2	500	2	1k	2	2k	2	5k	2	20k
2.5	40	2.5	160	2.5	400	2.5	800	2.5	1.6k	2.5	4k	2.5	16k
4	25	4	100	4	250	4	500	4	1k	4	2.5k	4	10k
5	20	5	80	5	200	5	400	5	800	5	2k	5	8k
10	10	10	40	10	100	10	200	10	400	10	1k	10	4k
20	5	20	20	20	50	20	100	20	200	20	500	20	2k
50	2	50	8	50	20	50	40	50	80	50	200	50	800
100	1	100	4	100	10	100	20	100	40	100	100	100	400

- Spindle     - Spindle Speed     - Factor    K - 1000

HA Series Viscometer													
1		2		3		4		5		6		7	
.5	400	.5	1.6k	.5	4k	.5	8k	.5	16k	.5	40k	.5	160k
1	200	1	800	1	2k	1	4k	1	8k	1	20k	1	80k
2	100	2	400	2	1k	2	2k	2	4k	2	10k	2	40k
2.5	80	2.5	320	2.5	800	2.5	1.6k	2.5	3.2k	2.5	8k	2.5	32k
5	40	5	160	5	400	5	800	5	1.6k	5	4k	5	16k
10	20	10	80	10	200	10	400	10	800	10	2k	10	8k
20	10	20	40	20	100	20	200	20	400	20	1k	20	4k
50	4	50	16	50	40	50	80	50	160	50	400	50	1.6k
100	2	100	8	100	20	100	40	100	80	100	200	100	800

- Spindle     - Spindle Speed     - Factor    K - 1000

HB Series Viscometer													
1		2		3		4		5		6		7	
.5	1.6k	.5	6.4k	.5	16k	.5	32k	.5	64k	.5	160k	.5	640k
1	800	1	3.2k	1	8k	1	16k	1	32k	1	80k	1	320k
2	400	2	1.6k	2	4k	2	8k	2	16k	2	40k	2	160k
2.5	320	2.5	1.28k	2.5	3.2k	2.5	6.4k	2.5	12.8k	2.5	32k	2.5	128k
5	160	5	640	5	1.6k	5	3.2k	5	6.4k	5	16k	5	64k
10	80	10	320	10	800	10	1.63k	10	3.2k	10	8k	10	32k
20	40	20	160	20	400	20	800	20	1.6k	20	4k	20	16k
50	16	50	64	50	160	50	320	50	640	50	1.6k	50	6.4k
100	8	100	32	100	80	100	160	100	320	100	800	100	3.2k

- Spindle     - Spindle Speed     - Factor    K - 1000

**Dial reading x Factor = Viscosity in cP (mPa•s)**

Example: LVT Viscometer with #1 spindle at 6 rpm

Dial Reading: 75    Factor: 10

$75 \times 10 = 750 \text{ cP (mPa}\cdot\text{s)}$

Full scale viscosity range for any speed and spindle combination is equal to the factor x 100.

**Factor x 100 = Full scale range**

Example: LVT Viscometer with #1 spindle at 6 RPM

Full Scale Range:  $10 \times 100 = 1,000 \text{ cP}$

### II.3 Considerations for Making Measurements

In taking viscosity measurements with the Dial Viscometer there are two considerations which pertain to the low viscosity limit of effective measurement.

- 1) Viscosity measurements should be accepted within the equivalent % Torque Range from 10% to 100% for any combination of spindle/speed rotation.
- 2) Viscosity measurements should be taken under laminar flow conditions, not under turbulent flow conditions.

The first consideration has to do with the precision of the instrument. All Dial Viscometers have a full scale range precision of  $\pm 1\%$  for any spindle/speed combination. We discourage taking readings below 10% of range because the potential viscosity error of  $\pm 1\%$  is a relatively high number compared to the instrument reading.

The second consideration involves the mechanics of fluid flow. All rheological measurements of fluid flow properties should be made under laminar flow conditions. Laminar flow is flow wherein all particle movement is in layers directed by the shearing force. For rotational systems, this means all fluid movement must be circumferential. When the inertial forces on the fluid become too great, the fluid can break into turbulent flow wherein the movement of fluid particles becomes random and the flow can not be analyzed with standard math models. This turbulence creates a falsely high viscometer reading with the degree of non-linear increase in reading being directly related to the degree of turbulence in the fluid.

For the following geometries, we have found that an approximate transition point to turbulent flow occurs as follows:

- 1) No. 1 LV Spindle: 15 cP at 60 RPM
- 2) No. 1 RV Spindle: 100 cP at 50 RPM
- 3) UL Adapter: 0.85 cP at 60 RPM

Turbulent conditions will exist in these situations whenever the RPM/cP ratio exceeds the values listed above.

## APPENDIX B - Viscosity Ranges

LV Series Viscometers with Spindles #1 - #4

RV/HA/HB Series Viscometers with Spindles #1 - #7

Viscosity Range (cP)		
Viscometer	Minimum	Maximum
LVF	15	100,000
LVT	15	2,000,000
RVF	100	2,000,000
RVT	100	8,000,000
HAT	200	16,000,000
HBT	800	64,000,000

# Rotoviskometer *Brookfield-Cone And Plate* Bagian 1

## APPENDIX A - Cone/Plate Viscometer Set-Up

The Wells-Brookfield Dial Viscometer uses the same operating instruction procedures as described in this manual, however, the gap between the cone and the plate must be mechanically adjusted before measurements are made. This is done by moving the plate (built into the sample cup, Part No. CP-44Y) up towards the cone until two small pins (one in the cone, the second mounted on the plate) contact slightly, and then by separating (lowering) the plate 0.0005inch (0.013mm).

Note that the Wells-BrookfieldCone/Plate Viscometer requires the use of a circulating temperature bath controlling temperature to within  $\pm 0.1$  °C. The bath is connected to the ports on the CP-44Y sample cup. Brookfield offers a complete line of refrigerated and non-refrigerated (tap water cooling) circulating water baths.

The following example assumes that the C/P Viscometer is set up on the Model A laboratory stand and has been leveled. Cone CP-40 is used in the illustrations. The water bath is on, set at the same test temperature that you will use when measuring the viscosity of your sample and is connected to the sample cup (Figure A1). Allow enough time for the bath to reach the test temperature. The Viscometer speed knob should be set at 10 or 12 rpm with motor OFF.


Figure A1

Remove the sample cup and attach the cone to the Viscometer (Note: left hand threads), using the spindle wrench to hold the Viscometer shaft (Figure A2). Note: Lift up gently on the spindle wrench.


Figure A2

Attach the sample cup and swing the clip under the cup to secure it in place. Take care to avoid hitting the cone with the cup (Figure A3). Note: Do not add test sample during the gap setting procedure.


Figure A3

Start the Viscometer running at 10 or 12 rpm. Readings of % Torque (dial reading) will be utilized for setting the cone/plate gap. If the display reading jumps to 0.5% of scale (or higher), or will not settle to zero (indicating that the pins on the cone and cup are hitting), unscrew the adjustment ring by turning it to the left (clockwise as you look down on the instrument) until the reading settles at 0.0% (Figure A4).


Figure A4

Turn the adjustment ring to the right (counterclockwise as you look down on the instrument) in small increments (one or two divisions of the ring) while watching the Viscometer dial (Figure A5). Remember that you are trying to raise the cup so that the pin in the cone touches the pin in the cup. Once you have found this hit point, you can "back off" to create the desired gap between cone and cup. *You must wait at least 6 seconds between movements of the ring.*

Turn the adjustment ring until the display reading jumps from 0.0 to 0.5% (HA and HB models) or from 0.0 to 0.5 % (RV models) or 0.0 to 1.0% (LV models). This is the HIT POINT. The cup may have to be adjusted up and down several times before this deflection is satisfactorily reached. The LV series Viscometer will be the *most sensitive*.


Figure A5

When you are satisfied the pins are just hitting (by observing acceptable % readings), make a pencil mark on the adjustment ring directly under the index mark on the pivot housing (Figure A6).


Figure A6

Turn the adjustment ring to the left (clockwise as you look down on the instrument) exactly the width of one division from the pencil mark you have just made. This will separate the pins by 0.0005" (0.013 mm). This is a very important step because if it is not done, the pins may touch constantly and may cause wear and damage to the cone, cup and instrument. The Viscometer is now mechanically set and ready for sample measurement (Figure A7).


Figure A7

Each of the five available cones has a specific sample volume, as shown in Table A1. Note that the correct amount of sample fluid should cover the cone face and up over the edge less than 1mm as shown in (Figure A8).


Figure A8

- Notes:
- The cup may be removed without resetting the gap.
  - Remove the spindle from the viscometer when you clean it.
  - Find the hit point every time the spindle is attached.

<u>Cone</u>	<u>Sample Volume</u>	<u>Cone Angle</u>
CP-40	0.5 ml	0.8°
CP-41	2.0 ml	3.0°
CP-42	1.0 ml	1.565°
CP-51	0.5 ml	1.565°
CP-52	0.5 ml	3.0°

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
# VISKOMETER VT-04F RION


## Specifications

	VT-03F	VT-04F
Measurement range	No. 4 rotor: 2 to 33 mPa·s No. 5 rotor: 15 to 150 mPa·s No. 3 rotor: 50 to 300 mPa·s	No. 3 rotor: 0.3 to 13 dPa·s (with No. 3 cup) No. 1 rotor: 3 to 150 dPa·s (with JIS 300 mL beaker*2) No. 2 rotor: 100 to 4000 dPa·s (with JIS 300 mL beaker*2)
Sample fluid capacity	Approx. 480 mL (with Cup A or Cup B)	No. 1 and No. 2 rotor (with JIS 300 mL beaker*2) approx. 300 mL No. 3 rotor (with No. 3 cup) approx. 170 mL Clearance between rotor end and cup bottom: about 15 mm
Measurement accuracy (and reproducibility for VT-04F)	Within ±5 % of scale maximum (using supplied cup A or cup B) Scale calibrated according to JIS Z 8809 : 2000 viscosity calibration standard	Measurement accuracy ±10 % of indicated value Reproducibility ±5 % (calibrated with silicon oil)
Rotor speed	62.5 rpm	
Power requirements	6 V DC (four IEC R6P batteries) or AC adapter VA-05A or B*1 Current consumption at maximum torque approx. 100 mA	6 V DC (four IEC R6P batteries) or AC adapter VA-05A or B*1 Current consumption at maximum torque approx. 300 mA
Dimensions	98 (W) × 181 (D) × 40 (H) mm (without protruding parts)	
Weight	Approx. 570 g (without batteries)	
Supplied accessories	No. 3 rotor (dia. 45 × 47 × 160 mm) SUS304 1 No. 4 rotor (dia. 78 × 46 × 159 mm) A1050 (alumite) 1 No. 5 rotor (dia. 61.2 × 36 × 149 mm) A1050 (alumite) 1 Cup A (dia. 92 × 76 mm) A1050 (alumite) 1 Cup B (dia. 92 × 76 mm) A1050 (alumite) 1 IEC R6P (size AA) battery 4 Instruction manual 1 Inspection certificate 1	No. 1 rotor (dia. 24 × 53 × 166 mm) SUS304 1 No. 2 rotor (dia. 15 × 1 × 113 mm) SUS304 1 No. 3 rotor (dia. 45 × 47 × 160 mm) SUS304 1 No. 3 Cup (dia. 52.6 × 75 mm) SUS304 1 Rotor extension (900 mm · 300× 3) SUS304 1 IEC R6P (size AA) battery 4 Instruction manual 1 Inspection certificate 1
Optional accessory	AC adapter VA-05A or B*, Stand VA-04	

\*1 VA-05A, 120 V VA-05B, 220 V \*2 JIS R 3503 : 1994, φ78×103 (H)

### Rotors and Cups (unit: mm)

\*1 Material: Stainless steel \*2 Material: Aluminate \*3 Cup B has a 30 mm diameter hole at the bottom. Use Cup B to measure viscosity by directly by dipping Rotor in the fluid.

### Sample amount for measurement

	VT-03F	VT-04F
Cup A	approx.460 mL	—
No.3 Cup	—	approx.170 mL
Commercially available 300 mL beaker	—	approx.350 mL

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
Pustaka:

Sinko PJ, 2012, *Physical Pharmacy and Pharmaceutical Sciences* 6<sup>th</sup> ed., Lippincott Williams and Wilkins. Page 415.

**Fig. 17-3.** Caking diagram, showing the flocculation of a bismuth subnitrate suspension by means of the flocculating agent monobasic potassium phosphate. (From A. Martin and J. Swarbrick, in *Sprows' American Pharmacy*, 6th Ed., Lippincott, Philadelphia, 1966, p. 205. With permission.)

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

Pustaka : Drug Stability for Pharmaceutical Scientist page 72-75.

### 3.1.4 $\beta$ -Lactam Antibiotics


$\beta$ -Lactam antibiotics are antibacterial agents that contain a  $\beta$ -lactam ring in their molecular structure, and they include penicillins, cephalosporins, monobactams, and carbapenems:


Penicillins

Cephalosporins


Monobactams


Carbapenems

The  $\beta$ -lactam ring is a four-member lactam ring that undergoes hydrolysis. In general, lactams are, like amides, chemically stable, but due to the structural strain of the four-member ring, the  $\beta$ -lactam ring is relatively rapidly hydrolyzed in aqueous solutions. The stability of the lactam ring is affected by the substituents (R). Hydrolysis of the  $\beta$ -lactam ring is catalyzed by intramolecular catalysis under acidic conditions:

Substituents that pull electrons away from the oxygen (i.e., electron withdrawing groups) decrease the rate of hydrolysis, while substituents that push electrons toward the oxygen (i.e., electron donating groups) enhance the rate. While ampicillin and amoxicillin are relatively stable under acidic conditions and can be given orally, benzylpenicillin is very unstable and is rapidly hydrolyzed in the stomach, which results in an extremely low oral bioavailability. *Amoxicillin* can be found in solid dosage forms (i.e., tablets and capsules), as oral suspensions, and as solutions for parenteral administration. Amoxicillin has four ionization forms and the following degradation pathway is based on the pH-rate profile (Fig. 3.10):


Figure 3.10 Sketch of the pH-rate profile for hydrolysis of amoxicillin in aqueous buffer solution at 35° C. The observed rate constant has the units of  $h^{-1}$  [12].

Based on Fig. 3.10, the equation for  $k_{obs}$  can be written as:

$$k_{obs} = k_H[H^+]f_{AH_3} + k'_H[H^+]f_{AH_2} + k'_0f_{AH_2} + k'_{OH}[OH^-]f_{AH_2} \\ + k''_{OH}[OH^-]f_{AH} + k'''_{OH}[OH^-]f_A$$

$$-\frac{d[A]_T}{dt} = k_{obs}[A]_T$$

$$[A]_T = [AH_3] + [AH_2] + [AH] + [A]$$

$$f_{AH_3} = \frac{[H^+]^3}{[H^+]^3 + [H^+]^2K_{a1} + [H^+]K_{a1}K_{a2} + K_{a1}K_{a2}K_{a3}}$$


$$f_{AH_2} = \frac{[H^+]^2K_{a1}}{[H^+]^3 + [H^+]^2K_{a1} + [H^+]K_{a1}K_{a2} + K_{a1}K_{a2}K_{a3}}$$

$$f_{AH} = \frac{[H^+]K_{a1}K_{a2}}{[H^+]^3 + [H^+]^2K_{a1} + [H^+]K_{a1}K_{a2} + K_{a1}K_{a2}K_{a3}}$$

$$f_A = \frac{K_{a1}K_{a2}K_{a3}}{[H^+]^3 + [H^+]^2K_{a1} + [H^+]K_{a1}K_{a2} + K_{a1}K_{a2}K_{a3}}$$

Maximum stability is at a pH of about 5, at which the shelf-life of an aqueous amoxicillin solution is less than one week at room temperature. Due to its instability, amoxicillin is marketed, not as ready-made solutions, but as a granular powder for preparation of aqueous suspensions for oral administration and as a lyophilized powder for reconstitution of a parenteral solution. Moisture in solid state formulations (e.g., tablets, capsules, and powders) promotes hydrolysis of the  $\beta$ -lactam ring.


Polyols (i.e., compounds with several hydroxyl groups like glucose, sucrose, dextrans, glycerol, mannitol, and sorbitol) catalyze  $\beta$ -lactam hydrolysis in aqueous solutions [13,14] through the formation of chemically unstable  $\beta$ -lactam/polyol complexes [15]. Thus,  $\beta$ -lactam antibiotics should not be dissolved in or mixed with carbohydrates or dissolved in, for example, isotonic glucose solutions. Likewise, metal ions such as  $\text{Cu}^{2+}$ ,  $\text{Zn}^{2+}$ , and  $\text{Mn}^{2+}$ , catalyze the cleavage of the  $\beta$ -lactam ring through formation of complexes [16],

where  $k_c$  is much larger than  $k_f$ . For example, for benzylpenicillin and  $\text{Cu}^{2+}$ ,  $K_c$  was determined to be  $190 \text{ M}^{-1}$ ,  $k_f$  to be  $0.15 \text{ M}^{-1} \text{ s}^{-1}$ , and  $k_c$   $1.2 \times 10^7 \text{ M}^{-1} \text{ s}^{-1}$  ( $k_c/k_f = 8 \times 10^7$ ) at  $30^\circ \text{C}$ . Thus, minute amounts of  $\text{Cu}^{2+}$  ions can have a significant destabilizing effect on  $\beta$ -lactam antibiotics.

### 3.1.5 Nitrogen Mustards

Nitrogen mustards are an important class of alkylating anticancer drugs that owe their activity to their ability to form reactive aziridinium ring by intramolecular displacement of the chloride by the amine nitrogen. In the presence of water, the reactive aziridinium ring (i.e., the ethylene-immonium ion) reacts with water to form inactive degradation products:

#### Pustaka :

Connors KA, Amidon GL, Kennon L, 1997, *Chemical Stability of Pharmaceuticals*, John Wiley & Sons. Page 129-150.